

Sonifier Cell Disruptor

**SFX 150 Microtips****OVERVIEW**

Microtips for the SFX Sonifier are available in four sizes to handle volumes up to 150 ml. The microtips are also designed to process samples in the smallest diameter vessels.

The 1/8" microtip is already provided with the SFX Sonifiers.

Part #	Processing Volume	Tip Diameter	Approx. Length (working)
109-122-1066	1ml & under	3/32" (2.4mm)	4"
109-122-1065	3 -10 ml	1/8" (3.2mm)	4"
109-122-1182	5 -25 ml	3/16" (4.8mm)	4"
109-122-1064	10 - 150 ml	1/4" (6.4mm)	4"
109-116-1566	3-10 ml (4) tips	1/8" (3.2mm)	4"

Note: The 109-116-1566 has (4) 1/8" microtips on a 1" center horn.

Note: These tips are for use with the Branson 150 watt SFX and SLP Sonifier only.

**Americas**

Branson Ultrasonics Corp.
41 Eagle Road
Danbury, CT 06813-1961
T: 203-796-0400
F: 203-796-9838
www.Emerson.com/Branson

Europe

Branson Ultraschall
Niederlassung der Emerson
Technologies GmbH & Co. OHG
Waldstrasse 53-55
63128 Dietzenbach, Germany
T: +49-6074-497-0
F: +49-6074-497-199
www.Emerson.com/Branson

Asia

Branson Ultrasonics (Shanghai) Co., Ltd.
758 Rong Le Dong Road
Song Jiang, Shanghai, PRC, 201613
T: 86-21-3781-0588
F: 86-21-5774-5100
www.Emerson.com/Branson